




VITABLOCS® TriLux® for CEREC®



the versatile, three-layer ceramic blocks for the CEREC® System



VITA



VITABLOCS® TriLux for CEREC® are manufactured from the VITABLOCS Mark II, the fine-structure feldspar ceramic that is kind to grinding tools and has been clinically proven 9 million times over, and demonstrate wear characteristics corresponding to those of natural tooth substance. Thanks to a special production procedure it has become possible – in addition to the outstanding light transmission effect and the whitish fluorescence of the Mark II ceramic – to integrate 3 different degrees of colour saturation (chroma), and therefore 3 different degrees of translucency in one block. This in turn makes it possible

- to simulate the colour effects characteristic of a natural tooth with regard to translucency and shade intensity,
- to achieve better integration of the restoration into the patient's remaining natural tooth substance or natural dentition as appropriate.

The product was developed specifically for dental use and must be processed according to the instructions for use. The information/illustrations contained herein do not constitute an assurance of properties and are not binding. Date of issue of these instructions for use: 03-2005.

VITABLOCS® TriLux for CEREC® are certified and bear the CE mark **CE 0124**

CEREC® is a registered trade mark of Sirona Dental Systems GmbH

The individual layers

- The middle layer (body) corresponds to a normal degree of translucency.
- The top layer (enamel) is the least shade intensive and at the same time more translucent.
- The bottom layer (neck) is the most highly pigmented and the least translucent.

Indications

VITABLOCS® TriLux for CEREC® are particularly suitable for manufacturing partial crowns, anterior and posterior crowns as well as veneers.

Assortments

VITABLOCS® TriLux for CEREC® are available in the sizes TRI-12 (10x12x15 mm) in packs of 10 pcs and TRI-14 (12x14x18 mm) in packs of 5 pcs. Both sizes are available in the three following VITA SYSTEM 3D-MASTER® shades for CEREC: 1M2C, 2M2C and 3M2C.

VITA

VITA Zahnfabrik H. Rauter GmbH & Co. KG
Postfach 1338 · D-79704 Bad Säckingen · Germany
Tel. +49/7761/562-222 · Fax +49/7761/562-446
www.vita-zahnfabrik.com · info@vita-zahnfabrik.com